


ABBE Refractometers

2 AWAJ Series



This model is a simple Abbe refractometer with thermometer and is extremely reliable. You can evaluate liquid samples, solid samples as well as paste samples. Key features of this refractometer are its sturdiness and its ease of use.

As an option there is also a beautifully-made aluminium transport and storage box available. The refraction index nD is measured.

- The models are universal, analogue Abbe refractometers
- The practical and robust design means that it can be used simply, efficiently and on a permanent basis, day in day out
- The integrated scale allows it to be used in many application areas and offers the best possible security for being able to read the measuring result accurately

Application :

- Main areas of application:
- Sugar industry (raw cane sugar)
- Pharmaceuticals
- Drinks industry
- Food processing industry
- Chemical industry
- Oil industry/Refineries
- Laboratories
- Training

Standard Accessories

- 1 x Abbe refractometer ABBE-REF 1
- 1 x screwdriver
- 1 x thermometer
- 1 x user manual

Model	2 AWAJ
Refractive index	1.300 ... 1.700 nD
Sugar content	0.0 ... 95 % Brix
Temperature	0 ... +70 °C
Resolution and accuracy	
Refractive index	0.0002 nD
Sugar content	0.1 % Brix
Scale division	
Refractive index	0.0005 nD
Sugar content	0.25 % Brix
Indicator/Split scale for dual measurement of	% Brix / nD
Dimension (L x W x H)	270 x 180 x 370mm
Weight/Approx.	3 kg

WORLDWIDEMARKETING

ACZET PVT. LTD.

Unit E2, Plot No. 15, WICEL Estate, Opp Seepz Gate Nol. 1, Andheri (E), Mumbai 400 093. (INDIA)
 Tel. : +91-22-4243 7700 • Fax : +91-22-4243 7800 • E-mail : sales@aczel.com • Web : www.aczel.com