

BIO SAFETY CABINET

MODEL: CLE-147



BIO-SAFETY CABINET

Ultra Clean Air Work Station, Vertical type, made of Pest & weather resistant board which is clad in laminated sheet Work stage thoroughly lined with Stainless Steel (S.S. - 304 or 316) with illumination by fluorescent tube lights. The unit is complete with H.E.P.A. Filtered clean exhaust duct. Efficiency up to 99.97%, retaining particles up to the size of 0.03 micron. To work on 220 230 V AC 50 Hz.

Available in following sizes:

- * 2'x2'x2'
- * 3'x2'x2'
- * 4'x2'x2'

Important Points for Quality

- | | |
|---|--|
| 1) Thickness of Plywood | : 19mm |
| 2) Thickness of Sunmyca | : 1mm |
| 3) Brand of H.E.P.A Filter | : Athlete (With Test Certificate) |
| 4) U.V. Tube | : 2 Nos. , (Brand: Phillips / Equivalent) |
| 5) Fluorescent tube | : Normal 40W Tube Light
(Brand: Philips / Equivalent) |
| 6) S. S. Top of Working Table | : 304 Grade |
| 7) Switches | : Heavels / Anchor / Bentex / L & T |
| 8) Electronic Motor | : ¼ H.P., ISI Marked, (Brand: BKP / Equivalent)- 1 No. |
| 9) Motor Protection Switch | : Isolation Switch Set of Good Quality
(Brand: Heavels / Anchor / Equivalent) |
| 10) Statics Pressure Manometer (Chemically Filled) Good Quality | : - 1 No. |