



**MODEL : HL-23**

Specification of **COSLAB** Model HL-23

- Viewing Head : Siedentopf Binocular/ Trinocular 30° Head, IPD 55-75 mm
- Eyepiece : Wide field Eyepiece WF 10X/ 18 mm
- Objectives : Anti Fungal Achromatic DIN Objectives 4X, 10X, 40X & 100X (Oil Immersion)
- Nosepiece : Quadruple Reverse Nosepiece (Turret)
- Mechanical Stage : Double Layers 140 X 140 mm/ 75 X 50 mm DOUBLE SLIDE HOLDER
- Condenser : NA1.25 Abbe Condenser with Iris Diaphragm & Filter
- Focusing System : Smooth Coaxial Coarse and fine Motion
- Illumination : 6V 20W Halogen Lamp or 3W LED, Brightness Adjustable.
- Dimension & weight : 36x26x46.5 cm & 8.5 kg.

**Optional Accessories :**

1. Fluorescent Attachment FLM-202
2. Phase Contrast Attachment PH-101 & Dark Field Attachment DF-101
3. Digital Camera or Digital Camera Eye piece with inbuilt USB Wire or CCD Camera.
4. **COSLAB** Image Analysis Software.
5. LED Illumination : Rechargeable Battery Back up more than 50 Hours.
6. Electronic Green Pointer Attachment.

Designs, Specifications and Picture are subject to change without any prior notice.

*Manufacturers & Exporters*

**COSMO LABORATORY EQUIPMENT**

3474, Shantipath (Rambagh Road), Ambala Cantt-133 001 (INDIA)

Tele/Fax : +91-171-2618278, 4006702 • Mb.: 098121-83773 / 094161-13230

E-mail : [info@coslabindia.com](mailto:info@coslabindia.com), [coslabindia@gmail.com](mailto:coslabindia@gmail.com) • Log on to : [www.coslabindia.com](http://www.coslabindia.com), [www.microscopesmanufacturer.com](http://www.microscopesmanufacturer.com)