

Basic Compact Scale CG Series



aczel



Standard Features

- Compact and space saving.
- Rugged design ABS plastic housing with efficient bench space usage.
- Simple to used.
- Auto External Calibration.
- Backlight LCD Display for easy viewing.
- Weighing Unit gm, dwt, oz, ozt, available
- Auto zero tracking
- AC or DC Powersuppl : Outfit AC / DC adopter.
- Over Load Indication.
- Improved repeatability for better results.
- **GSM Unit Available only CG 302MG Series.**

Compact Scale for basic weighing applications

Model	CG 302MG	CG 302M	CG 502M	CG 1001M	CG 3001M	CG 5001M
Capacity	300gm	300 gm	500 gm	1000 gm	3000 gm	5000 gm
Accuracy	0.01 gm	0.01 gm	0.01 gm	0.1 gm	0.1 gm	1 gm
Pan Size (mm/inch)	118 Ø / 4.7"Ø			155 x 155 mm / 6.1" x 6.1"	120 x 120 mm / 4.7" x 4.7"	155 x 155 mm / 6.1" x 6.1"
Weighing Units	GSM	Gram, oz, dwt, ozt				
Tare Range	Full					
Calibration	Auto					
Stability Time	1.5 Second					
Auto Power OFF	Yes					
Temperature	10 - 45°C					
Display	LCD with Backlight					

Specification are subject to change without prior notice.
Images shown are for illustrations purpose only.

**WORLDWIDEMARKETING
ACZET PVT. LTD.**

Unit E2, Plot No. 15, WICEL Estate, Opp Seepz Gate No. 1, Andheri (E), Mumbai 400 093. (INDIA)
Tel. : +91-22-4243 7700 • Fax : +91-22-4243 7800 • E-mail : sales@aczel.com • Web : www.aczel.com