


Inclined Binocular Microscope



MODEL: HL-9

- Viewing Head : Binocular / Trinocular 45° rotatable 360°
- Eyeiece : Wide Field Eyeiece 10X
- Objectives : Achromatic objectives 4X, 10X, 40X & 100X
- Nosepiece : Quadruple Nosepiece
- Mechanical Stage : Coaxial Mechanical Stage.
- Focusing Range : Separate coarse and fine adjustment
- Illumination : 6V 20W Halogen Lamp, Brightness Adjustable.

Optional Accessories:

1. Digital Camera or Digital Camera Eye piece with inbuilt USB Wire or CCD Camera.
2. COSLAB Image Analysis Software.