# PRECISION BALANCES CY SERIES





CITIZEN's new and improved precision series incorporates unique and unequalled High end technology built to match the highest international standards. The result is TEN fold increase in weighing performance which maintains our proven, simple and easy to use operator interface. Attractive quality of this new design is demonstrated by the CE mark of conformity along with the CITIZEN's ISO 9001: 2000 certification.







Citizen precision balances come with all the intelligent capability you need for your daily routine application other than "Basic Weighing" like



- Piece Counting
- Percentage (%) weighing
- Formulation
- Totalisation
- · CY-C Series precision balances comes with fully Automatic Internal Calibration.



• Density measurement for solid and liquid [with optional eqp & software]

11:22AM

CG 3102 9821161 r0.1.04 1234567

7654321

11:2301

50 pcs 1999.97 g 9999.85 ct

50 pcs

......

PRINT FORMATS

22-Mar-10

Calibration:

Model

N-ID

Diff.

22-Mar-10

Name:

Citizen

External Cal Done Diff. 0.00s

11:13AM

CG 3182

9821161 r8.1.04 1234567

External

+ 2000.009

11:13AM

- · Fully Automatic Internal calibration with various 'PSAC' Following features. Enables user for highly perfect and Precise weighing
- INTIAL Calibration
- · Self Calibration according to change in the ambient temperature

22-Mar-18

Model Ser.no. Ver.no. ID

LID:

nRef uRef

22-Mar-18

- User selected TIME Calibration
- Formulation
- Pipette Calibration
- Below weighing
- Animal Weighing



























### STANDARD FEATURES

High contrast, Large Backlite LCD display for easy viewing with A.E.P. (Advanced Eye Protection)

Standard RS 232 Cinterface/

Hanger for Below Balance weighing Automatic External calibration

Compiles GLP/GMP Standard.

Dye cast aluminum design for long term stability and

accurate results

Various weighing units like gm, mg, ct, oz, dwt, mom, GN,g, lb, ozt, t1T, t1H, t1S, Bat, MS.

User selectable stability and filter level

USB Interface available (Optional)

#### Statistical Calculation Function:

The balance display and / or output statistical calculation data as No. (number of data), SUM (total), MAX (maximum), MIN (minimum), R [range (maximur

num - minimum)j, Av E (average), SD (standard deviation) and CV (coefficient of variation).								
data:	Standard Level		Professional Le	vel				

Technical data	tandard Level			Professional Level								
MODEL	CY 120	CY 220	CY 320	CY 420	CY 510	CY 720	CY 1003	CY 220C	CY 320C	CY 510C	CY 1003C	
CAPACITY	120 gm	220 gm	320 gm	420 gm	510 gm	720 gm	1000 gm	220 gm	320 gm	510 gm	1000 gm	
READABILITY	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	
REPEATABILITY (sd)	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.001 gm	0.002 gm	0.001 gm	0.001 gm	0.001 gm	0.002 gm	
LINEARITY (+/-)	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	0.002 gm	
PAN SIZE (mm/inch) 128 x 128 / 5.0" x 5.0"												
Response Time 2-3 sec												
DISPLAY	Back Lite LCD Display											
CALIBRATION				Automatic External			Motorised Internal					
UNITS OF MEASURE G.mg,ct,GN,mo,oz,dwt,t1T,t1H,t1S,t1S,mom,Bat,MS												
TARE RANGE	ARE RANGE Full											
OPERATING TEMP.	OPERATING TEMP. 15° C to 40° C											
HOUSING DIM.( WxDxH)		342mm x 212mm x 27					x 271mm				210 x 300 x 145mm 4.1 kg / 9.5 lb	
WEIGHT APPROX 5.2kg/ 10.5lb												
POWER SUPPLY	A/C Adapter 100 V-300 V +/- 20% 50-60 Hz											
Additional equipment:												
PRINTER CPR 02 UNDER WEIGHING HOOK												
DENSITY DETERMINATION KIT WEIGHTRACK COMPUTER SOFTWARE												
PC KEYBOARD CALIBRATION WEIGHT E2, F1, F2, M1 ( NABL, NPL, RRSL, STANDARD )								))				
ADDITIONAL LCD DISPLA		VIBROSORB-ANTI VIBRATION TABLE										
			V	VODI DW	IDEMARI	CETING:						

## WORLDWIDEMARKETING:

# CITIZEN SCALE (I) PVT. LTD.

Unit No. E2, Plot No. 11, WICEL Estate, Opp. Seepz Gate No.1, Andheri (E), Mumbai 400 093. (INDIA) Tel.: +91-22-4243 7700 • Fax: +91-22-4243 7800 • E:mail: sales@citizenscales.com • Web: www.citizenscales.com